

# Perifix<sup>®</sup> FX Springwound Epidural Catheter

Designed to enhance patient comfort and minimize risks of complications



Easy placement, optimized drug spread, and rare incidence of paresthesia and vascular placement.

# B. Braun Perifix® FX Springwound Epidural Catheter



With an inner **spring embedded** in polyamide nylon, the Perifix FX catheter provides the ideal balance of softness and stability.

**The result?** A design intended to provide safer, more accurate insertions that benefit patients and anesthesiologists alike. And it's available in a preferred multi-port design to ensure optimal drug spread, even under low-flow conditions.

### Perifix FX Advantages:

- Reduction of paresthesia and vessel puncture rates
- Flexibility
- Ease of placement and removal
- Consistency
- Reliability

**The bottom line?** Our advanced springwound epidural catheters are designed to enable easier insertion, enhanced patient comfort and minimized risk of paresthesia and vessel puncture.

READ WHAT USERS THINK ABOUT THE PERIFIX FX AT [www.perifixfx.bbraunusa.com/FX](http://www.perifixfx.bbraunusa.com/FX)



Flashback windows (distal and proximal)



Easy insertion and removal with reduced threat of breakage or uncoiling with the spring embedded in polyamide



Easier to advance to minimize the incidence of catheter repositioning for better analgesia quality




Centimeter markings


Transition from rigid body to flexible tip



Closed-tip, multi-port design with three lateral holes for optimal drug spread, fewer unilateral blocks and a potentially higher detection rate of intravascular catheter placement



Perifix FX is radiopaque



Also available in an open-tip, uni-port design

Reduced paresthesia and vessel puncture with a soft and atraumatic distal tip



Consistent feel and performance with minimal variation in tip deflection



**“Threads smoothly, rare vascularization of catheter, rare migration of catheter to subdural space. So far it's the best, most functional and easily placed catheter I've used.”**

**–G. Guffey, CRNA**

For more information and product trials, and to read what other users think about the Perifix® FX, see:  
[www.perifixfx.bbraunusa.com/FX](http://www.perifixfx.bbraunusa.com/FX)

**B | BRAUN**  
SHARING EXPERTISE

[www.bbraunusa.com](http://www.bbraunusa.com)

PC3-08 5/10 KFD 2M